

Bsa I GMP-grade (20 U/ μ L)

Product description

This product is a type IIs restriction endonuclease derived from the recombinant protein encoded by the BsaI gene in *Bacillus sphaericus* expressed by *E.coli*. Its recognition sequence is 5'-GGTCTCN1/N5-3'. Use to digest plasmids to prepare poly(A/T/G/C)-terminated linearized DNA fragments to obtain specific cohesive ends.

This product is produced in accordance with GMP process requirements and provided in a liquid form.

Specifications

Expression Host	Recombinant <i>E. coli</i> with Bsa I gene
Reaction Temperature	37°C
Storage Buffer	10mM Tris-HCl, 0.2M NaCl, 0.1mM EDTA, 1mM DTT, 50% Glycerol, 0.2mg/ml OsrHSA pH 7.4 \pm 0.2 (25°C)
Unit Definition	1 unit: The amount of enzyme required to digest 1 μ g of substrate DNA within 1 h at 37°C in a 50 μ L system.
Application	1. Digest the plasmid to prepare a linearized DNA fragment at the end of Poly (A/T/C/G); 2. Digestion of DNA to obtain specific sticky ends; 3. Linearize plasmid template before in-vitro transcription.

Components

Components No.	Name	10661ES03 (1 KU)	10661ES10 (10 KU)	10661ES60 (100 KU)
10661	Bsa I GMP-grade (20 U/ μ L)	50 μ L	500 μ L	5 mL

Storage

This product should be stored at -25 ~ -15°C for two years.

Instructions

Experimental methods

50 μ L reaction system

This step is suitable for linearization of 1 μ g DNA (\geq 100 nt) and can be scaled up according to experimental needs.

1. Add the following components in sequence:

Components	Volume
Plasmid DNA	1-2 μ g
10 \times Digestion Buffer 4	5.0 μ L
Bsa I (20 U/ μ L)	1.0 μ L

RNase-free ddH₂O

Up to 50 μ L

【Note】 10 \times Digestion Buffer 4(Cat#10668ES): 500 mM Potassium Acetate, 200 mM Tris-acetate, 100 mM Magnesium Acetate, 1 mg/ml OsrHSA, pH7.9@25°C

2. Incubate at 37°C 1 h;
3. DNA linearization is complete, and subsequent experiments can be performed.

Notes

1. Heat inactivation condition: incubate at 80°C for 20min.
2. Please operate with lab coats and disposable gloves, for your safety.