

## BsaI GMP-grade (20 U/ $\mu$ L)

### Product description

This product is a type IIS restriction endonuclease derived from the recombinant protein encoded by the BsaI gene in *Bacillus sphaericus* expressed by *E. coli*. Its recognition sequence is 5'-GGTCTCN1/N5-3'. Use to digest plasmids to prepare poly(A/T/G/C)-terminated linearized DNA fragments to obtain specific cohesive ends.

This product is produced in accordance with GMP process requirements and provided in a liquid form.

### Specifications

|                      |   |
|----------------------|---|
| Expression Host      | Recombinant <i>E. coli</i> with BsaI gene   |
| Reaction Temperature | 37°C  |
| Storage Buffer       | 10 mM Tris-HCl, 0.2 M NaCl, 0.1 mM EDTA, 1mM DTT, 50% Glycerol  |
| Unit Definition      | 1 unit: The amount of enzyme required to digest 1 $\mu$ g of substrate DNA within 1 h at 37°C in a 50 $\mu$ L system.   |
| Application          | 1. Digest the plasmid to prepare a linearized DNA fragment at the end of Poly (A/T/G/C);<br>2. Digestion of DNA to obtain specific sticky ends;<br>3. Linearize plasmid template before in-vitro transcription. |

### Components

| Components No. | Name                           | 10661ES03<br>(1 KU) | 10661ES10<br>(10 KU) | 10661ES60<br>(100 KU) |
|----------------|--------------------------------|---------------------|----------------------|-----------------------|
| 10661          | BsaI GMP-grade (20 U/ $\mu$ L) | 50 $\mu$ L          | 500 $\mu$ L          | 5 mL                  |

### Storage

This product should be stored at -25 ~ -15°C for two years.

### Instructions

Experimental methods

50  $\mu$ L reaction system

This step is suitable for linearization of 1  $\mu$ g DNA ( $\geq$ 100 nt) and can be scaled up according to experimental needs.

1. Add the following components in sequence:

| Components            | Volume      |
|-----------------------|-------------|
| Plasmid DNA           | 1-2 $\mu$ g |
| 10×Digestion Buffer 4 | 5.0 $\mu$ L |
| BsaI (20 U/ $\mu$ L)  | 1.0 $\mu$ L |

**【Note】** 10× Digestion Buffer 4 (Cat#10668ES): 500 mM Potassium Acetate, 200 mM Tris-acetate, 100 mM Magnesium Acetate, 1 mg/ml OsrHSA, pH7.9@25°C

2. Incubate at 37° C 1 h;
3. DNA linearization is complete, and subsequent experiments can be performed.

### Notes

1. Heat inactivation condition: incubate at 80° C for 20min.
2. Please operate with lab coats and disposable gloves, for your safety.